


WOUND CARE PRODUCTS

Advanced Wound Management

Calcium Alginate Dressings



Algicell Calcium Alginate are soft and white sterile dressings that absorb moderate to large amounts of wound drainage and may help control minor bleeding. Alginate ropes can be used to fill or pack the dead space in a wound. A secondary dressing is needed to hold alginates in place and to protect the wound from outside contaminants.

CATALOG NUMBER	Description	Units/Case
88022	Algicell Calcium Alginate Dressing- Sterile 2" x 2"	40 units/case
88044	Algicell Calcium ALginate Dressing -Sterile 4" x 4"	40 units/case
88048	Algicell Calcium Alginate Dressing- Sterile 4" x 8"	40 units/case
88112	Algicell Alginate Rope Dressing- Sterile, 3/4" x 12"	40 units/case
88022	Algicell Calcium Alginate Dressing- Sterile 2" x 2"	40 units/case
88044	Algicell Calcium ALginate Dressing -Sterile 4" x 4"	40 units/case
88048	Algicell Calcium Alginate Dressing- Sterile 4" x 8"	40 units/case
88112	Algicell Alginate Rope Dressing- Sterile, 3/4" x 12"	40 units/case