


WOUND CARE PRODUCTS

Advanced Wound Management

Hydrocolloid Dressing



PrimaCol Hydrocolloid Dressing mix with wound exudate to form a soft mass that allows for easy removal while protecting the wound from outside contamination from bacteria, fecal matter or urine.

CATALOG NUMBER	Description	Units/Case
85144	PrimaCol Hydrocolloid Dressing 4" x 4"	80 units/case
85166	PrimaCol Hydrocolloid Dressing 6" x 6"	80 units/case
85188	PrimaCol Hydrocolloid Dressing 8" x 8"	80 units/case