






*Science is at the core of our business.*

Our commitment to excellence has made Victus one of the world's most respected names in nutritional and healthcare products.

Healthcare professionals and health-conscious people, young and old alike, look to our leadership for quality institutional and consumer products.



## Medical Products - Dosi-Flow I.V. Administration and Extension Sets

| Catalog Number | Description  | Image  | Case Quantity |
|----------------|--|--|---------------|
| 259-0250       | Dosi-Flow 1 Non-DEHP Two-piece female luer lock connector, flow range of 5 to 250 ml/h, one latex free injection site 6" (15.24 cm) from the distal end and two-piece spin male luer lock connector. Length 18" (45.68 cm), priming volume 2.5 ml.   |    | 50 per case   |
| 259-0270       | Dosi-Flow 1 Non-DEHP Two-piece female luer lock connector, flow range of 5 to 250 ml/h, one Victus needle free injection site 6" (15.24 cm) from the distal end and two-piece spin male luer lock connector. Length 18" (45.68 cm), priming volume 2.5 ml.   |    | 50 per case   |
| 259-0253       | Dosi-Flow 3 Universal Adult Administration Non- DEHP Set with Standard Injection Site Macrodropper (20 drops/ml) universal vented set, macrodrip chamber with 15 micron filter, flow range of 5 to 250 ml/h and one latex free injection site 6" (15.24 cm) from the distal end. Length 83" (210.65 cm), priming volume 14 ml.                             |  | 50 per case   |
| 259-0293       | Dosi-Flow 3 Universal Adult Administration Non- DEHP Set with Victus Injection Site, Macrodropper (20 drops/ml) universal vented set, macrodrip chamber with 15 micron filter, flow range of 5 to 250 ml/h and one Victus needle free injection site 6" (15.24 cm) from the distal end. Length 83" (210.65 cm), priming volume 14 ml.                      |  | 50 per case   |
| 259-0300       | Dosi-Flow 3 Universal Adult Administration Non-DEHP Set with Victus Injection Site, Macrodropper (20 drops/ml) universal vented set, macrodrip chamber with 15 micron filter, flow range of 5 to 250 ml/h and one Victus needle free injection site 6" (15.24 cm) from the distal end, 0.22 micron filter. Length 84" (213.36 cm), priming volume 14.2 ml. |  | 50 per case   |